

Plate Star P

Application and Features

Buffered Developer Powder for pre-coated offset plates. Customers with fast plates should also take **Developer Inhibitor liquid**, Art.-No. 106244. The **Plate Star P** works self-cleaning with a high stability for a very very easy production with a Mixing System like **Variomix II**.

Usage

Temperature: 22 – 24⁰C
 Speed: minimum 20 sec.
 Dwelling Time.

Replenishment

Spray Bath Processors: 100 - 150 ml/m²
 Dive Bath Processors: 60 - 100 ml/m²

Packaging

Plate Star P

1 Package for 10 x 20 litres Developer

Art.-No. 106240

Developer Inhibitor liquid

1 Package with 10 x 100 ml Bottles

Art.-No.106244

Variomix II

Mixing System with 2 x 40 litres Tanks.
 For Developer and Fixing Powders.

Art.-No. 101383